



GGT5 (heavy chain, Cleaved-Gly387) mouse mAb

Catalog No	YP-mAb-04383
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GGT5 GGTLA1
Immunogen	Synthesized peptide derived from human GGT5 (heavy chain, Cleaved-Gly387)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat GGT5 (heavy chain, Cleaved-Gly387, protein was cleaved amino acid sequence between 387-388)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GGT5;GGT5 GGTLA1;GGT5 (heavy chain;Cleaved-Gly387);Gamma-glutamyltransferase 5 (GGT 5;EC 2.3.2.2;Gamma-glutamyl transpeptidase-related enzyme;GGT-rel;Gamma-glutamyltransferase-like activity 1;Gamma-glutamyltranspeptidase 5;Glutathione hydrolase 5;EC 3.4.19.13;Leukotriene-C4 hydrolase;EC 3.4.19.14) [Cleaved into: Gamma-glutamyltransferase 5 heavy chain;Gamma-glutamyltransferase 5 light chain]; Cleaved-Gly387)
Observed Band	42 65kD
Calculated Molecular Weight	42 65kD
Cell Pathway	Membrane ; Single-pass type II membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images