



ITIH3 (Cleaved-Asp651) mouse mAb

Catalog No	YP-mAb-04385
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ITIH3
Immunogen	Synthesized peptide derived from human ITIH3 (Cleaved-Asp651)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat ITIH3 (Cleaved-Asp651, protein was cleaved amino acid sequence between 651-652)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ITIH3;ITIH3 (Cleaved-Asp651);Inter-alpha-trypsin inhibitor heavy chain H3 (ITI heavy chain H3;ITI-HC3;Inter-alpha-inhibitor heavy chain 3;Serum-derived hyaluronan-associated protein;SHAP)
Observed Band	68 95kD
Calculated Molecular Weight	68 95kD
Cell Pathway	Secreted.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images