



KLK1/2/3 (Cleaved-Ile25) mouse mAb

Catalog No	YP-mAb-04387
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KLK1/2/3
Immunogen	Synthesized peptide derived from human KLK1/2/3 (Cleaved-Ile25)
Specificity	This antibody detects endogenous levels of Human KLK1/2/3 (Cleaved-Ile25, protein was cleaved amino acid sequence between 24-25)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KLK1/2/3;KLK1/2/3 (Cleaved-Ile25);Kallikrein-2 (EC 3.4.21.35;Glandular kallikrein-1;hGK-1;Tissue kallikrein-2)
Observed Band	25kD
Calculated Molecular Weight	25kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images