



MMP26 (Cleaved-Thr90) mouse mAb

Catalog No	YP-mAb-04393
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MMP26
Immunogen	Synthesized peptide derived from human MMP26 (Cleaved-Thr90)
Specificity	This antibody detects endogenous levels of Human MMP26 (Cleaved-Thr90, protein was cleaved amino acid sequence between 89-90)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MMP26;MMP26 (Cleaved-Thr90);Matrix metalloproteinase-26 (MMP-26;EC 3.4.24.-;Endometase;Matrilysin-2)
Observed Band	17 25kD
Calculated Molecular Weight	17 25kD
Cell Pathway	Secreted, extracellular space, extracellular matrix.
Tissue Specificity	Expressed specifically in uterus and placenta. Is also widely expressed in malignant tumors from different sources as well as in diverse tumor cell lines.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images