



NOTCH3 (Cleaved-Asp1628) mouse mAb

Catalog No	YP-mAb-04396
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NOTCH3
Immunogen	Synthesized peptide derived from human NOTCH3 (Cleaved-Asp1628)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat NOTCH3 (Cleaved-Asp1628, protein was cleaved amino acid sequence between 1628-1629)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NOTCH3;NOTCH3 (Cleaved-Asp1628);Neurogenic locus notch homolog protein 3 (Notch 3) [Cleaved into: Notch 3 extracellular truncation;Notch 3 intracellular domain]
Observed Band	78 260kD
Calculated Molecular Weight	78 260kD
Cell Pathway	Cell membrane ; Single-pass type I membrane protein.; [Notch 3 intracellular domain]: Nucleus. Following proteolytical processing NICD is translocated to the nucleus.
Tissue Specificity	Ubiquitously expressed in fetal and adult tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images