



PRSS8 (heavy chain, Cleaved-Ile45) mouse mAb

Catalog No	YP-mAb-04398
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRSS8
Immunogen	Synthesized peptide derived from human PRSS8 (heavy chain, Cleaved-Ile45)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat PRSS8 (heavy chain, Cleaved-Ile45, protein was cleaved amino acid sequence between 44-45)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRSS8;PRSS8 (heavy chain;Cleaved-Ile45);Prostasin (EC 3.4.21.-;Channel-activating protease 1;CAP1;Serine protease 8) [Cleaved into: Prostasin light chain;Prostasin heavy chain]; Cleaved-Ile45)
Observed Band	32kD
Calculated Molecular Weight	32kD
Cell Pathway	[Prostasin]: Cell membrane; Single-pass membrane protein.; [Prostasin light chain]: Secreted, extracellular space. Found in the seminal fluid. Secreted after cleavage of its C-terminus.; [Prostasin heavy chain]: Secreted, extracellular space. Found in the seminal fluid. Secreted after cleavage of its C-terminus.
Tissue Specificity	Found in prostate, liver, salivary gland, kidney, lung, pancreas, colon, bronchus and renal proximal tubular cells. In the prostate gland it may be synthesized in epithelial cells, secreted into the ducts, and excreted into the seminal fluid.
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images