



DYR1A Monoclonal Antibody

Catalog No	YP-mAb-04891
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DYRK1A DYRK MNB MNBH
Immunogen	Synthesized peptide derived from human protein . at AA range: 270-350
Specificity	DYR1A Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DYRK1A DYRK MNB MNBH;Dual specificity tyrosine-phosphorylation-regulated kinase 1A (EC 2.7.12.1) (Dual specificity YAK1-related kinase) (HP86) (Protein kinase minibrain homolog) (MNBH) (hMNB)
Observed Band	83
Calculated Molecular Weight	
Cell Pathway	Nucleus . Nucleus speckle .
Tissue Specificity	Ubiquitous. Highest levels in skeletal muscle, testis, fetal lung and fetal kidney.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images