



USF1 Monoclonal Antibody

Catalog No	YP-mAb-04923
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	USF1 BHLHB11 USF
Immunogen	Synthesized peptide derived from human protein . at AA range: 90-170
Specificity	USF1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	USF1 BHLHB11 USF;Upstream stimulatory factor 1 (Class B basic helix-loop-helix protein 11) (bHLHb11) (Major late transcription factor 1)
Observed Band	34
Calculated Molecular Weight	
Cell Pathway	Nucleus.
Tissue Specificity	Kidney,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images