



UBP28 Monoclonal Antibody

Catalog No	YP-mAb-04924
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	USP28 KIAA1515
Immunogen	Synthesized peptide derived from human protein . at AA range: 650-730
Specificity	UBP28 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	USP28 KIAA1515;Ubiquitin carboxyl-terminal hydrolase 28 (EC 3.4.19.12) (Deubiquitinating enzyme 28) (Ubiquitin thioesterase 28) (Ubiquitin-specific-processing protease 28)
Observed Band	118
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleoplasm .
Tissue Specificity	Brain, Eye,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images