



AP2A Monoclonal Antibody

Catalog No	YP-mAb-04947
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TFAP2A AP2TF TFAP2
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	AP2A Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TFAP2A AP2TF TFAP2;Transcription factor AP-2-alpha (AP2-alpha) (AP-2 transcription factor) (Activating enhancer-binding protein 2-alpha) (Activator protein 2) (AP-2)
Observed Band	48
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Prostate,Teratocarcinoma,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images