



EGR2 Monoclonal Antibody

Catalog No	YP-mAb-04965
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EGR2 KROX20
Immunogen	Synthesized peptide derived from human protein . at AA range: 340-420
Specificity	EGR2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AT591;CMT1D;CMT4E;E3 SUMO-protein ligase EGR2;E3 SUMO-protein transferase ERG2;EGR2;EGR2 KROX20;E3 SUMO-protein ligase EGR2 (EC 6.3.2.-) (AT591) (Early growth response protein 2) (EGR-2) (Zinc finger protein Krox-20)
Observed Band	60-70 kDa
Calculated Molecular Weight	476 aa, 50 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Endometrium,
Function	disease:Defects in EGR2 are a cause of Charcot-Marie-Tooth disease type 1D (CMT1D) [MIM:607678]. CMT1D is a form of Charcot-Marie-Tooth disease, the most common inherited disorder of the peripheral nervous system. Charcot-Marie-Tooth disease is classified in two main groups on the basis of electrophysiologic properties and histopathology: primary peripheral demyelinating neuropathy or CMT1, and primary peripheral axonal neuropathy or CMT2. Neuropathies of the CMT1 group are characterized by severely reduced nerve conduction velocities (less than 38 m/sec), segmental demyelination and remyelination with onion bulb formations on nerve biopsy, slowly progressive distal muscle atrophy and weakness, absent deep tendon



reflexes, and hollow feet.,disease:Defects in EGR2 are a cause of congenital hypomyelination neuropathy (CHN) [MIM:605253]. Inheritance can be autosomal dominant or recessive. R

Background

The protein encoded by this gene is a transcription factor with three tandem C2H2-type zinc fingers. Defects in this gene are associated with Charcot-Marie-Tooth disease type 1D (CMT1D), Charcot-Marie-Tooth disease type 4E (CMT4E), and with Dejerine-Sottas syndrome (DSS). Multiple transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Oct 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images