



ZN496 Monoclonal Antibody

Catalog No	YP-mAb-05015
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF496 ZKSCAN17
Immunogen	Synthesized peptide derived from human protein . at AA range: 10-90
Specificity	ZN496 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NIZP1;zinc finger protein 496;ZKSCAN17;ZNF496;ZN496;ZNF496 ZKSCAN17; Zinc finger protein 496 (Zinc finger protein with KRAB and SCAN domains 17)
Observed Band	67 kDa
Calculated Molecular Weight	587 aa, 67 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Colon,Lung,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 1 SCAN box domain.,similarity:Contains 5 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 1 SCAN box domain.,similarity:Contains 5 C2H2-type zinc fingers.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images