



TRI35 Monoclonal Antibody

Catalog No	YP-mAb-05033
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRIM35 HLS5 KIAA1098
Immunogen	Synthesized peptide derived from human protein . at AA range: 320-400
Specificity	TRI35 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	E3 ubiquitin-protein ligase TRIM35;EC:2.3.2.27;HLS5;KIAA1098;tripartite motif containing 35;TRI35;TRIM35 HLS5 KIAA1098;Tripartite motif-containing protein 35 (Hemopoietic lineage switch protein 5)
Observed Band	45~50 kDa
Calculated Molecular Weight	57 kDa
Cell Pathway	Cytoplasm . Nucleus . Found predominantly in cytoplasm with a granular distribution. Found in punctuate nuclear bodies (By similarity). .
Tissue Specificity	Brain,Liver,Lung,
Function	domain:The RING finger domain and the coiled-coil region are required for the apoptosis-inducing activity.,function:May play a role as a tumor suppressor. Implicated in the cell death mechanism.,similarity:Belongs to the TRIM/RBCC family.,similarity:Contains 1 B box-type zinc finger.,similarity:Contains 1 B30.2/SPRY domain.,similarity:Contains 1 RING-type zinc finger.,subcellular location:Found predominantly in cytoplasm with a granular distribution. Found in punctuate nuclear bodies.,

**Background**

The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The function of this protein has not been identified. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images