



COX1 Monoclonal Antibody

Catalog No	YP-mAb-05051
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MT-CO1 COI COXI MTCO1
Immunogen	Synthesized peptide derived from human protein . at AA range: 380-460
Specificity	COX1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COX1;COI;COX 1;COX I;COXI;MT-CO1 COI COXI MTCO1;Cytochrome c oxidase subunit 1 (EC 1.9.3.1) (Cytochrome c oxidase polypeptide I)
Observed Band	56kD
Calculated Molecular Weight	56kD
Cell Pathway	Mitochondrion inner membrane ; Multi-pass membrane protein .
Tissue Specificity	Blood,Bone fossil,Bones,Breast cancer,Distant normal tissue,Glioma,Para-can
Function	catalytic activity:4 ferrocytochrome c + O(2) + 4 H(+) = 4 ferricytochrome c + 2 H(2)O.,disease:Defects in MT-CO1 are a cause of anemia sideroblastic acquired idiopathic (AISA) [MIM:516030]; a disease characterized by inadequate formation of heme and excessive accumulation of iron in mitochondria., disease:Defects in MT-CO1 are a cause of cytochrome c oxidase deficiency (COX deficiency) [MIM:220110]; also called mitochondrial complex IV deficiency. COX deficiency is a clinically heterogeneous disorder. The clinical features are ranging from isolated myopathy to severe multisystem disease, with onset from infancy to adulthood.,disease:Defects in MT-CO1 are a cause of Leber hereditary optic neuropathy (LHON) [MIM:535000]. LHON is a maternally inherited disease resulting in acute or subacute loss of central vision, due to optic nerve dysfunction. Cardiac conduction defects and neurological d



Background

catalytic activity: $4 \text{ ferrocytochrome } c + \text{O}(2) + 4 \text{ H}(+) = 4 \text{ ferricytochrome } c + 2 \text{ H}(2)\text{O}$.,disease:Defects in MT-CO1 are a cause of anemia sideroblastic acquired idiopathic (AISA) [MIM:516030]; a disease characterized by inadequate formation of heme and excessive accumulation of iron in mitochondria.,disease:Defects in MT-CO1 are a cause of cytochrome c oxidase deficiency (COX deficiency) [MIM:220110]; also called mitochondrial complex IV deficiency. COX deficiency is a clinically heterogeneous disorder. The clinical features are ranging from isolated myopathy to severe multisystem disease, with onset from infancy to adulthood.,disease:Defects in MT-CO1 are a cause of Leber hereditary optic neuropathy (LHON) [MIM:535000]. LHON is a maternally inherited disease resulting in acute or subacute loss of central vision, due to optic nerve dysfunction. Cardiac conduction defects and neurological defects have also been described in some patients. LHON results from primary mitochondrial DNA mutations affecting the respiratory chain complexes.,disease:Defects in MT-CO1 are associated with recurrent myoglobinuria [MIM:550500]. Myoglobinuria consists of excretion of myoglobin in the urine.,function:Cytochrome c oxidase is the component of the respiratory chain that catalyzes the reduction of oxygen to water. Subunits 1-3 form the functional core of the enzyme complex. CO I is the catalytic subunit of the enzyme. Electrons originating in cytochrome c are transferred via the copper A center of subunit 2 and heme A of subunit 1 to the bimetallic center formed by heme A3 and copper B.,pathway:Energy metabolism; oxidative phosphorylation.,similarity:Belongs to the heme-copper respiratory oxidase family.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images