



KRT37 Monoclonal Antibody

Catalog No	YP-mAb-05122
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KRT37 HHA7 HKA7 KRTHA7
Immunogen	Synthesized peptide derived from human protein . at AA range: 310-390
Specificity	KRT37 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KRT37;KRT37 HHA7 HKA7 KRTHA7;Keratin;type I cuticular Ha7 (Hair keratin; type I Ha7) (Keratin-37) (K37); type I cuticular Ha7 (Hair keratin; type I Ha7) (Keratin-37) (K37)
Observed Band	49kD
Calculated Molecular Weight	49kD
Cell Pathway	intermediate filament,extracellular exosome,
Tissue Specificity	Scalp,
Function	miscellaneous:There are two types of hair/microfibrillar keratin, I (acidic) and II (neutral to basic).,similarity:Belongs to the intermediate filament family.,
Background	The protein encoded by this gene is a member of the keratin gene family. As a type I hair keratin, it is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images