



HCN3 Monoclonal Antibody

Catalog No	YP-mAb-05962
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HCN3 KIAA1535
Immunogen	Synthesized peptide derived from human protein . at AA range: 150-230
Specificity	HCN3 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HCN3;KIAA1535;HCN3 KIAA1535;Potassium/sodium hyperpolarization-activated cyclic nucleotide-gated channel 3
Observed Band	86 kDa
Calculated Molecular Weight	774 aa 86 kDa
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Detected in brain.
Function	domain:The segment S4 is probably the voltage-sensor and is characterized by a series of positively charged amino acids at every third position., function:Putative hyperpolarization-activated ion channel exhibiting weak selectivity for potassium over sodium ions.,similarity:Belongs to the potassium channel HCN family.,similarity:Contains 1 cyclic nucleotide-binding domain., subunit:The potassium channel is probably composed of a homo- or heterotetrameric complex of pore-forming subunits.,
Background	This gene encodes a multi-pass membrane protein that functions as a voltage gated cation channel. The encoded protein is a member of a family of closely related cyclic adenosine monophosphate-binding channel proteins. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images