



ABC3A Monoclonal Antibody

Catalog No	YP-mAb-05974
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	APOBEC3A
Immunogen	Synthesized peptide derived from human protein . at AA range: 80-160
Specificity	ABC3A Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	APOBEC3A;Probable DNA dC->dU-editing enzyme APOBEC-3A (EC 3.5.4.-) (Phorbolin-1)
Observed Band	21
Calculated Molecular Weight	
Cell Pathway	Nucleus. Cytoplasm.
Tissue Specificity	Expressed in peripheral leukocytes with higher expression in CD14-positive phagocytic cells. Highly expressed in keratinocytes and in periphery blood monocytes. Also detected in non-lymphoid tissues including lung and adipose tissues. Found at high levels in colorectal adenocarcinoma, Burkitt's lymphoma and chronic myelogenous leukemia.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images