



SAAL1 Monoclonal Antibody

Catalog No	YP-mAb-05993
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SAAL1
Immunogen	Synthesized peptide derived from human protein . at AA range: 320-400
Specificity	SAAL1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein SAAL1;serum amyloid A like 1;SPACIA1;Synoviocyte proliferation-associated in collagen-induced arthritis protein 1;SAAL1
Observed Band	54-55 kDa
Calculated Molecular Weight	475 aa, 54 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Highly expressed in testis and ovary, and to a lesser extent in the lung, spleen and the heart (at protein level).
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the SAAL1 family.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the SAAL1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images