



RFTN2 Monoclonal Antibody

Catalog No	YP-mAb-06011
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RFTN2 C2orf11
Immunogen	Synthesized peptide derived from human protein . at AA range: 370-450
Specificity	RFTN2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C2orf11;Raft linking protein 2;Raftlin 2;raftlin family member 2;RFTN2;RFTN2 C2orf11;Raftlin-2 (Raft-linking protein 2)
Observed Band	56-60 kDa
Calculated Molecular Weight	501 aa, 56 kDa
Cell Pathway	Cell membrane ; Lipid-anchor .
Tissue Specificity	Brain,
Function	similarity:Belongs to the raftlin family.,
Background	similarity:Belongs to the raftlin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images