



RGF1C Monoclonal Antibody

Catalog No	YP-mAb-06031
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RASGEF1C
Immunogen	Synthesized peptide derived from human protein . at AA range: 100-180
Specificity	RGF1C Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ35841;OTTHUMP161522;Ras GEF domain containing family member 1C; RASGEF 1C;RasGEF domain containing family member 1C;RasGEF domain family member 1C;RGF1C;RASGEF1C;Ras-GEF domain-containing family member 1C
Observed Band	51kD
Calculated Molecular Weight	53 kDa; 53kD(Calculated).
Cell Pathway	intracellular,
Tissue Specificity	Brain,Testis,
Function	function:Guanine nucleotide exchange factor (GEF).,similarity:Contains 1 N-terminal Ras-GEF domain.,similarity:Contains 1 Ras-GEF domain.,
Background	function:Guanine nucleotide exchange factor (GEF).,similarity:Contains 1 N-terminal Ras-GEF domain.,similarity:Contains 1 Ras-GEF domain.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images