



SHE Monoclonal Antibody

Catalog No	YP-mAb-06162
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SHE
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SHE Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RP11 350G8.8;SHE;SH2 domain-containing adapter protein E
Observed Band	55 kDa
Calculated Molecular Weight	54 kDa
Cell Pathway	
Tissue Specificity	Spleen,
Function	similarity:Contains 1 SH2 domain.,
Background	similarity:Contains 1 SH2 domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images