



TMC7 Monoclonal Antibody

Catalog No	YP-mAb-06319
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMC7
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TMC7 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C2364B5;C6324K23Rik;DKFZp781O2274;FLJ2124;RGD1564267;TMC 7;TMC7;TMC7_HUMAN;Transmembrane channel like;Transmembrane channel like protein 7;Transmembrane channel-like protein 7
Observed Band	79kD
Calculated Molecular Weight	84 kDa; 84kD(Calculated).
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Brain,Colon,
Function	function:May function as ion channels, transporters, or modulators of such., similarity:Belongs to the TMC family.,
Background	function:May function as ion channels, transporters, or modulators of such., similarity:Belongs to the TMC family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images