



UBP35 Monoclonal Antibody

Catalog No	YP-mAb-06343
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	USP35 KIAA1372 USP34
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	UBP35 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.4.19.12;KIAA1372;Ubiquitin carboxyl-terminal hydrolase 35;Ubiquitin-specific-processing protease 35;USP34;UBP35;USP35 KIAA1372 USP34; Ubiquitin carboxyl-terminal hydrolase 35 (EC 3.4.19.12) (Deubiquitinating enzyme 35) (Ubiquitin thioesterase 35) (Ubiquitin-specific-processing protease 35)
Observed Band	84 kDa
Calculated Molecular Weight	1018 aa, 113 kDa
Cell Pathway	
Tissue Specificity	Expressed in testis, pancreas and skeletal muscle.
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., similarity:Belongs to the peptidase C19 family.,tissue specificity:Expressed in testis, pancreas and skeletal muscle.,
Background	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., similarity:Belongs to the peptidase C19 family.,tissue specificity:Expressed in testis, pancreas and skeletal muscle.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images