



ZBTB2 Monoclonal Antibody

Catalog No	YP-mAb-06404
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZBTB2 KIAA1483 ZNF437
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ZBTB2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1483;Zinc finger and BTB domain-containing protein 2;ZNF437;ZBTB2; ZBTB2 KIAA1483 ZNF437
Observed Band	56kD
Calculated Molecular Weight	7-75 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Melanoma,Stomach,
Function	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 4 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 4 C2H2-type zinc fingers.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images