



S39A6 Monoclonal Antibody

Catalog No	YP-mAb-06418
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC39A6 LIV1 ZIP6
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	S39A6 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LIV 1;LIV-1/ZIP6;SLC39A6;ZIP6;Estrogen-regulated protein LIV-1;S39A6;SLC39A6 LIV1 ZIP6;Zinc transporter ZIP6 (Estrogen-regulated protein LIV-1) (Solute carrier family 39 member 6) (Zrt- and Irt-like protein 6) (ZIP-6)
Observed Band	85-103 kDa, 48 and 32 kDa
Calculated Molecular Weight	85 kDa
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Highly expressed in the breast, prostate, placenta, kidney, pituitary and corpus callosum. Weakly expressed in heart and intestine. Also highly expressed in cells derived from an adenocarcinoma of the cervix and lung carcinoma.
Function	function:May act as a zinc-influx transporter.,induction:Up-regulated by estrogen in breast cancer cells lines.,PTM:N-glycosylated.,similarity:Belongs to the ZIP transporter (TC 2.A.5) family.,tissue specificity:Highly expressed in the breast, prostate, placenta, kidney, pituitary and corpus callosum. Weakly expressed in heart and intestine. Also highly expressed in cells derived from an adenocarcinoma of the cervix and lung carcinoma.,

**Background**

Zinc is an essential cofactor for hundreds of enzymes. It is involved in protein, nucleic acid, carbohydrate, and lipid metabolism, as well as in the control of gene transcription, growth, development, and differentiation. SLC39A6 belongs to a subfamily of proteins that show structural characteristics of zinc transporters (Taylor and Nicholson, 2003 [PubMed 12659941]).[supplied by OMIM, Mar 2008]

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images