



ALKB5 Monoclonal Antibody

Catalog No	YP-mAb-06679
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ALKBH5 ABH5 OFOXD1
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ALKB5 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ALKBH5 ABH5 OFOXD1;Probable alpha-ketoglutarate-dependent dioxygenase ABH5 (EC 1.14.11.-) (Alkylated DNA repair protein alkB homolog 5)
Observed Band	43
Calculated Molecular Weight	
Cell Pathway	Nucleus speckle .
Tissue Specificity	Widely expressed, with highest expression in lung, followed by testis, pancreas, spleen and ovary.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images