



H2B1B Monoclonal Antibody

Catalog No	YP-mAb-06760
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HIST1H2BB H2BFF
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	H2B1B Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HIST1H2BB H2BFF;Histone H2B type 1-B (Histone H2B.1) (Histone H2B.f) (H2B/f)
Observed Band	13
Calculated Molecular Weight	
Cell Pathway	Nucleus. Chromosome.
Tissue Specificity	Epithelium,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images