



HPLN1 Monoclonal Antibody

Catalog No	YP-mAb-06930
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HAPLN1 CRTL1
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	HPLN1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HAPLN1 CRTL1;Hyaluronan and proteoglycan link protein 1 (Cartilage-linking protein 1) (Cartilage-link protein) (Proteoglycan link protein)
Observed Band	38
Calculated Molecular Weight	
Cell Pathway	Secreted, extracellular space, extracellular matrix.
Tissue Specificity	Widely expressed. Weakly expressed in the brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images