



VWA3A Monoclonal Antibody

Catalog No	YP-mAb-07149
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VWA3A
Immunogen	Synthesized peptide derived from human protein . at AA range: 710-790
Specificity	VWA3A Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VWA3A;von Willebrand factor A domain-containing protein 3A
Observed Band	130kD
Calculated Molecular Weight	130kD
Cell Pathway	Secreted .
Tissue Specificity	Brain,Trachea,Uterus,
Function	similarity:Contains 2 VWFA domains.,
Background	similarity:Contains 2 VWFA domains.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images