



DEF1 Monoclonal Antibody

Catalog No	YP-mAb-07358
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DEFA1 DEF1 DEFA2 MRS; DEFA1B
Immunogen	Synthesized peptide derived from human protein . at AA range: 41-90
Specificity	DEF1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DEFA1 DEF1 DEFA2 MRS;DEFA1B;Neutrophil defensin 1 (Defensin; alpha 1) (HNP-1) (HP-1) (HP1) [Cleaved into: HP 1-56;Neutrophil defensin 2 (HNP-2) (HP-2) (HP2)]
Observed Band	10
Calculated Molecular Weight	
Cell Pathway	Secreted.
Tissue Specificity	Liver,Lung,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images