



CYFP1 Monoclonal Antibody

Catalog No	YP-mAb-07676
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CYFIP1 KIAA0068
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	CYFP1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CYFIP1 KIAA0068;Cytoplasmic FMR1-interacting protein 1 (Specifically Rac1-associated protein 1) (Sra-1) (p140sra-1)
Observed Band	137
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Cytoplasm, perinuclear region . Cell projection, lamellipodium . Cell projection, ruffle . Cell junction, synapse, synaptosome . Highly expressed in the perinuclear region (By similarity). Enriched in synaptosomes (By similarity). Also enriched in membrane ruffles and at the tips of lamellipodia (By similarity). .
Tissue Specificity	B-cell lymphoma,Bone marrow,Brain,Hepatoma,Lung,Placenta,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images