



PLF4 Monoclonal Antibody

Catalog No	YP-mAb-07711
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PF4 CXCL4 SCYB4
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	PLF4 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PF4 CXCL4 SCYB4;Platelet factor 4 (PF-4) (C-X-C motif chemokine 4) (Iroplact) (Oncostatin-A) [Cleaved into: Platelet factor 4; short form]
Observed Band	11
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Leukocyte,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images