



DFFB Monoclonal Antibody

Catalog No	YP-mAb-07736
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DFFB CAD DFF2 DFF40
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	DFFB Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DFFB CAD DFF2 DFF40;DNA fragmentation factor subunit beta (EC 3.-.-.) (Caspase-activated deoxyribonuclease) (CAD) (Caspase-activated DNase) (Caspase-activated nuclease) (CPAN) (DNA fragmentation factor 40 kDa subunit) (DFF-40)
Observed Band	37
Calculated Molecular Weight	
Cell Pathway	Cytoplasm. Nucleus.
Tissue Specificity	Brain,Fetal brain,Lung,Pancreas,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images