



IRK6 Monoclonal Antibody

Catalog No	YP-mAb-07748
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KCNJ6 GIRK2 KATP2 KCNJ7
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	IRK6 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KCNJ6 GIRK2 KATP2 KCNJ7;G protein-activated inward rectifier potassium channel 2 (GIRK-2) (BIR1) (Inward rectifier K(+) channel Kir3.2) (KATP-2) (Potassium channel; inwardly rectifying subfamily J member 6)
Observed Band	46
Calculated Molecular Weight	
Cell Pathway	Membrane; Multi-pass membrane protein.
Tissue Specificity	Most abundant in cerebellum, and to a lesser degree in islets and exocrine pancreas.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images