



DIAP2 Monoclonal Antibody

Catalog No	YP-mAb-07773
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DIAPH2 DIA
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	DIAP2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DIAPH2 DIA;Protein diaphanous homolog 2 (Diaphanous-related formin-2) (DRF2)
Observed Band	121
Calculated Molecular Weight	
Cell Pathway	[Isoform 3]: Cytoplasm, cytosol. Early endosome. Isoform 3 is cytosolic but when coexpressed with RHOD, the 2 proteins colocalize to early endosomes.
Tissue Specificity	Expressed in testis, ovary, small intestine, prostate, lung, liver, kidney and leukocytes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images