



DJC24 mouse mAb

Catalog No	YP-mAb-09188
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DNAJC24 DPH4 ZCSL3
Immunogen	Synthesized peptide derived from human DJC24 AA range: 39-89
Specificity	This antibody detects endogenous levels of human DJC24
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	173A21Rik;26127M2Rik;AV66965;AW24712;CSL type zinc finger containing protein 3;CSL-type zinc finger-containing protein 3;DJC24_HUMAN;DnaJ (Hsp4) homolog;subfamily C;member 24;DnaJ homolog subfamily C member 24;DNAJC24;DPH4 homolog (JJJ3;S. cerevisiae);DPH4 homolog;DPH4;JJJ3 homolog;JJJ3;MmDjC7;RGD156471;RP23-232D1.6;ZCSL3;Zinc finger;CSL domain containing 3;CSL type containing 3;DJC24;DNAJC24 DPH4 ZCSL3;DnaJ homolog subfamily C member 24 (CSL-type zinc finger-containing protein 3) (DPH4 homolog)
Observed Band	
Calculated Molecular Weight	17kD
Cell Pathway	Cytoplasm, cytoskeleton .
Tissue Specificity	
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images