



RFOX3 mouse mAb

Catalog No	YP-mAb-09226
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RBFOX3
Immunogen	Synthesized peptide derived from human RFOX3 AA range: 32-82
Specificity	This antibody detects endogenous levels of human RFOX3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ56884;FLJ58356;Fox-1 homolog C;fox1 homolog C;Fox3;FOX3NeuN; hexaribonucleotide binding protein 3;HRNBP3;NEUN;neuronal nuclei;Rbfox3; RFOX3_HUMAN;RNA binding protein fox-1 homolog 3;RNA binding protein;fox 1 homolog (C. elegans) 3;RFOX3;RNA binding protein fox-1 homolog 3 (Fox-1 homolog C)
Observed Band	
Calculated Molecular Weight	45-55kD; 34kD
Cell Pathway	Nucleus . Cytoplasm . Largely restricted to neuronal nuclei. However, significant cytoplasmic localization in neurons from brains from HIV-infected individuals with cognitive impairment. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images