



TCP1L mouse mAb

Catalog No	YP-mAb-09246
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TCP10L PRED77
Immunogen	Synthesized peptide derived from human TCP1L AA range: 68-118
Specificity	This antibody detects endogenous levels of human TCP1L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TCP1L;TCP10L PRED77;T-complex protein 10A homolog 2 (T-complex protein 10A-2) (TCP10A-2) (TCP10-like)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in liver and testis. Expressed in the seminiferous tubules (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images