



TTC6 mouse mAb

Catalog No	YP-mAb-09266
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TTC6
Immunogen	Synthesized peptide derived from human TTC6 AA range: 246-296
Specificity	This antibody detects endogenous levels of human TTC6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	213O14Rik;CGI 135;CGI135;FIS 1;FIS1;Fis1 homolog;S. cerevisiae;homolog of; FIS1_HUMAN;Fission 1 (mitochondrial outer membrane) homolog (S. cerevisiae);Fission 1 (mitochondrial outer membrane) homolog (yeast);Fission 1 (mitochondrial outer membrane) homolog;Fission 1 homolog;H NH132A1.6;hFis 1;hFis1;Mitochondrial fission 1 protein;mitochondrial fission molecule; Tetratricopeptide repeat domain 11;Tetratricopeptide repeat protein 11;TPR repeat protein 11;TTC 11;TTC6;Tetratricopeptide repeat protein 6 (TPR repeat protein 6)
Observed Band	
Calculated Molecular Weight	17 kDa; 17kD
Cell Pathway	
Tissue Specificity	
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images