



PCDB9 mouse mAb

Catalog No	YP-mAb-09337
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PCDHB9 PCDH3H
Immunogen	Synthesized peptide derived from human PCDB9 AA range: 443-493
Specificity	This antibody detects endogenous levels of human PCDB9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PCDB9;PCDHB9 PCDH3H;Protocadherin beta-9 (PCDH-beta-9) (Protocadherin-3H)
Observed Band	
Calculated Molecular Weight	88kD
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

