



ENT1 mouse mAb

Catalog No	YP-mAb-09379
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Synthesized peptide derived from human ENT1 AA range: 245-295
Specificity	This antibody detects endogenous levels of human ENT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ENT-1;Equilibrative NBMPR-sensitive nucleoside transporter;Equilibrative nitrobenzylmercaptapurine riboside-sensitive nucleoside transporter;es nucleoside transporter;Nucleoside transporter;es-type;ENT1;HERV-T_19q13.11 provirus ancestral Env polyprotein (Envelope polyprotein) (HERV-T Env protein) [Includes: Surface protein (SU);Transmembrane protein (TM)]
Observed Band	70 kDa
Calculated Molecular Weight	50KD
Cell Pathway	Cell membrane ; Single-pass membrane protein . At the origin, this retroviral envelope protein was localized in the virion. .
Tissue Specificity	Expressed at higher level in thyroid. Expressed at lower level in adrenal, bone marrow, brain, breast, kidney, ovary, placenta, prostate, skin, testis and trachea.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images