



BROMI mouse mAb

Catalog No	YP-mAb-09438
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BROMI C6orf170 C6orf171
Immunogen	Synthesized peptide derived from human BROMI AA range: 202-252
Specificity	This antibody detects endogenous levels of human BROMI
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACN;BANF;Bilateral acoustic neuroma;MERL_HUMAN;Merlin;Moesin ezrin radixin like protein;Moesin ezrin radixin like;Moesin-ezrin-radixin-like protein;Neurofibromatosis 2;Neurofibromatosis type 2;Neurofibromatosis2;Neurofibromin 2;Neurofibromin-2;Neurofibromin2;NF 2;Nf2;SCH;Schwannomerlin;Schwannomin;BROMI;BROMI C6orf170 C6orf171;Protein broad-minded
Observed Band	
Calculated Molecular Weight	7kDa; 7kD
Cell Pathway	Cytoplasm . Cell projection, cilium .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images