



AGAL mouse mAb

Catalog No	YP-mAb-09486
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GLA
Immunogen	Synthesized peptide derived from human AGAL AA range: 138-188
Specificity	This antibody detects endogenous levels of human AGAL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	beta galactosidase 2;BETA GALACTOSIDASE PROTECTIVE PROTEIN;beta-galactosidase 2;beta-galactosidase protective protein;betagalactosidase 2; Carboxypeptidase C;Carboxypeptidase L;carboxypeptidase Y-like kininase; Cathepsin A;Ctsa;deamidase;EC 3.4.16.5;Glactosialidosis;GLB2;Goldberg Syndrome;GSL;lysosomal carboxypeptidase A;Lysosomal protective protein 2 kDa chain;Lysosomal protective protein;Lysosomal protective protein deficiency; NEURAMINIDASE BETA GALACTOSIDASE EXPRESSION;NGBE; Neuraminidase deficiency with beta-galactosidase deficiency;OTTHUMP31778; OTTHUMP31781;PPCA;PPCA deficiency;PPGB;PPGB_HUMAN;Protective protein cathepsin A;Protective protein for beta galactosidase;Protective protein for beta-galactosidase;Protective protein/cathepsin A deficiency;urinary kininase; AGAL;GLA;Alpha-galactosidase A (EC 3.2.1.22) (Alpha-D-galactosidase A) (Alpha-D-galactoside galactohydrolase) (Melibiase) (Agalsidase)
Observed Band	
Calculated Molecular Weight	54kDa; 54kD
Cell Pathway	Lysosome.
Tissue Specificity	

**Function****Background****matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images