



GGTA1 mouse mAb

Catalog No	YP-mAb-09499
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GGTA1P GGTA1
Immunogen	Synthesized peptide derived from human GGTA1 AA range: 18-68
Specificity	This antibody detects endogenous levels of human GGTA1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CHM;Choroideraemia protein;Choroideremia;DXS54;FLJ38564;GGTA;HSD 32; MGC1271;Rab escort protein 1;Rab geranylgeranyltransferase component A; Rab proteins geranylgeranyltransferase component A 1;RAE1_HUMAN;REP 1; REP-1;REP1;TCD;TCD protein;GGTA1;GGTA1P GGTA1;Inactive N-acetyllactosaminide alpha-1;3-galactosyltransferase (Glycoprotein alpha-galactosyltransferase 1 pseudogene)
Observed Band	
Calculated Molecular Weight	73kD
Cell Pathway	Golgi apparatus, Golgi stack membrane ; Single-pass type II membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images