



SEN34 mouse mAb

Catalog No	YP-mAb-09505
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TSEN34 LENG5 SEN34
Immunogen	Synthesized peptide derived from human SEN34 AA range: 238-288
Specificity	This antibody detects endogenous levels of human SEN34
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:4.6.1.16;HsSen34;LENG5;PCH2C;SEN34;TSEN34 LENG5 SEN34;tRNA-splicing endonuclease subunit Sen34 (EC 3.1.27.9) (Leukocyte receptor cluster member 5) (tRNA-intron endonuclease Sen34) (HsSen34)
Observed Band	34-36 kDa
Calculated Molecular Weight	34 kDa
Cell Pathway	Nucleus . Nucleus, nucleolus . May be transiently localized in the nucleolus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images