



LAPM5 mouse mAb

Catalog No	YP-mAb-09569
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LAPTM5 KIAA0085
Immunogen	Synthesized peptide derived from human LAPM5 AA range: 179-229
Specificity	This antibody detects endogenous levels of human LAPM5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LAPM5_HUMAN;LAPTM5;Lysosomal-associated multitransmembrane protein 5; Lysosomal-associated transmembrane protein 5;Retinoic acid-inducible E3 protein;LAPM5;LAPTM5 KIAA0085;Lysosomal-associated transmembrane protein 5 (Lysosomal-associated multitransmembrane protein 5) (Retinoic acid-inducible E3 protein)
Observed Band	
Calculated Molecular Weight	15,25kD; 3kD
Cell Pathway	Lysosome membrane ; Multi-pass membrane protein .
Tissue Specificity	Preferentially expressed in adult hematopoietic tissues. High levels in lymphoid and myeloid tissues. Highly expressed in peripheral blood leukocytes, thymus, spleen and lung, followed by placenta, liver and kidney.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images