



RS16 mouse mAb

Catalog No	YP-mAb-09584
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPS16
Immunogen	Synthesized peptide derived from human RS16 AA range: 39-89
Specificity	This antibody detects endogenous levels of human RS16
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AUTSX 3;AUTSX3;DKFZp686A2416;Mbd 5;Mbd5;MECP 2;MeCP 2 protein; MeCP-2 protein;Mecp2;MECP2_HUMAN;Methyl CpG binding protein 2 (Rett syndrome);Methyl CpG binding protein 2;Methyl-CpG-binding protein 2;MRX 16; MRX 79;MRX16;MRX79;MRXS 13;MRXS13;MRXSL;PPMX;RS;RTS;RTT;WBP 1;WBP1;RS16;RPS16;40S ribosomal protein S16
Observed Band	
Calculated Molecular Weight	16kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

