



CD92 Polyclonal Antibody

Catalog No	YP-Ab-13917
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	SLC44A1
Immunogen	The antiserum was produced against synthesized peptide derived from human SLC44A1. AA range:391-440
Specificity	CD92 Polyclonal Antibody detects endogenous levels of CD92 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CD92;CDw92;Choline transporter-like protein 1;CHTL1;CTL1;SLC44A1;Solute carrier family 44 member 1;CD antigen CD92
Observed Band	70-73 kDa
Calculated Molecular Weight	73 kDa
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Mitochondrion outer membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in various cells of the hematopoietic system.
Function	function:Probable choline transporter. May be involved in membrane synthesis and myelin production.,similarity:Belongs to the CTL (choline transporter-like) family.,tissue specificity:Expressed in various cells of the hematopoietic system.,
Background	function:Probable choline transporter. May be involved in membrane synthesis and myelin production.,similarity:Belongs to the CTL (choline transporter-like) family.,tissue specificity:Expressed in various cells of the hematopoietic system.,
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images