




TCR β Polyclonal Antibody


Catalog No	YP-Ab-13964
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	TRBC1
Immunogen	The antiserum was produced against synthesized peptide derived from human TCR beta. AA range:66-115
Specificity	TCR β Polyclonal Antibody detects endogenous levels of TCR β protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	TCR β ;TRBC1;T-cell receptor beta-1 chain C region
Observed Band	20kD
Calculated Molecular Weight	20kD
Cell Pathway	Cell membrane .
Tissue Specificity	Lymph node: intermediate region,
Function	similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,
Background	similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TCR β Polyclonal Antibody



Western blot analysis of lysate from COLO205 cells treated with Forskolin, using TCR β antibody.