



K1210 mouse mAb

Catalog No	YP-mAb-09642
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KIAA1210
Immunogen	Synthesized peptide derived from human K1210 AA range: 747-797
Specificity	This antibody detects endogenous levels of human K1210
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DKFZp434K121;HD gene regulatory region binding protein 2;HD gene regulatory region-binding protein 2;HD-regulating factor 2;HDBP-2;HDBP2;HDRF-2;Huntington disease gene regulatory region-binding protein 2;Huntington's disease gene regulatory region binding protein 2;Papillomavirus regulatory factor 1;papillomavirus regulatory factor PRF 1;Papillomavirus-binding factor;PBF;PRF-1;PRF1;PTTG1IP;Si-1-8-14;Zinc finger protein 395;ZN395_HUMAN;znf395;K1210;KIAA1210;Uncharacterized protein KIAA1210
Observed Band	
Calculated Molecular Weight	55kD
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, acrosome . Localizes to the sex body in spermatocyte, acrosome, and near the ectoplasmic specialization. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

